

JAMMU AND KASHMIR MEDICAL SUPPLIES CORPORATION LTD.

(Public Sector Undertaking of the Government of Jammu and Kashmir)
Corporate Head Office: Plot No. 58, Friends Colony Satyam Road Trikuta Nagar Jammu
Corporate Office: Opp J&K Motor Garage Deptt near Hajj House Bemina Srinagar Telefax: 01942432008 emall: jkmscl@amail.com; jkmscl@amail.com</a

CORRIGENDUM

Subject: E Bid For the Finalization of Annual Rate Contract for the Procurement of "Suture Material" Uploaded vide No. JKMSCL/SUTURE Material/450 dated 03-03-2021

In continuation to e-bid uploaded for the finalization of Annual Rate contract for the procurement of "Suture Material" Vide No. JKMSCL/SUTURE Material/450 dated 03-03-2021 and pre-bid meeting held thereof on 17-03-2021, following amendments, in light of suggestions made by the various venders /Firms/Bidders, duly approved by the technical experts are hereby made as under:-

Item		Description			
Code	Sutures Materials	Needle type & size	Size	Length	Amendments
SM- 212	Braided silk suture	3/8 Circle Cutting Needle 60 MM	1-0	76cm (±10cm)	Please read No-1 in place of 1-0 and ±15cm in place of ±10cm
SM- 216	Braided silk suture	1/2 Circle Round Body/ Taper point 25- 26 MM	3-0	45cm (± 5cm)	Please read 90 in place of 45 cm
SM- 231	Braided silk suture	Pre-cut sutures Black Braided silk sterile	1	30cm (±5cm)	Please read 75cm in place of 30 cm
SM- 255	Polyglactin/Glycoli de &lactide braided absorbable suture (dyed)	5/8 Circle Cutting 16- 17 MM	4-0	75cm (±10cm)	Please read ±15cm in place of ±10cm
SM- 269	Polyglactin/Glycolid e & lactide braided Absorbable suture Undyed	3/8 Circle Cutting 26 MM	2-0	90cm (±10cm)	Please read ±15cm in place of ±10cm
SM- 282	Polyglactin/Glycolid e &lactide Braided triclosan Coated suture	1/2 circle reverse cutting OS 40 mm	2	90cm (±10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 283	Polyglactin/Glycolid e &lactide triclosan Braided Antibacterial suture (dyed)	1/2 Circle Round Body/Taper point17- 20 MM	4-0	70cm (±10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 284	Polyglactin /Glycolide &lactide, triclosan Braided Antibacterial suture (dyed)	1/2 Circle Round Body/Taper point (Heavy) 40 MM	1	90cm (± 10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 285	Polyglactin/Glycolid e &lactide, triclosan Braided Antibacterial suture (dyed)	1/2 Circle Round Body/Taper point 20- 22 MM	3-0	70cm (± 10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 286	Polyglactin/Glycolid e &lactide, triclosan Braided Antibacterial suture (dyed)	1/2 Circle Reverse Cutting 40 MM	1	90cm (± 10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 294	Nylon polyamide Monofilament	3/8 Circle Cutting Needle 22 MM	5-0	90cm (±10cm)	Please read Cutting Needle 12-22 MM in place of





JAMMU AND KASHMIR MEDICAL SUPPLIES CORPORATION LTD.

(Public Sector Undertaking of the Government of Jammu and Kashmair)
Corporate Head Office: Plot No. 58, Friends Colony Satyam Road Trikuta Nagar Jammu
Corporate Office: Opp J&K Motor Garage Deptt near Hajj House Bemina Srinagar Telefax: 01942432008 email: jkmscij@gmail.com; epm.jkmsci@gmail.com website: www.jkmscl.nic.in

CORRIGENDUM

Subject: E Bid For the Finalization of Annual Rate Contract for the Procurement of "Suture Material" Uploaded vide No. JKMSCL/SUTURE Material/450 dated 03-03-2021

"Suture Material" Vide No. JKMSCL/SUTURE Material/450 dated 03-03-2021 and pre-bid meeting held thereof on 17-03-2021, following amendments, in light of suggestions made by the various denders /Firms/Bidders, duly approved by the technical experts are hereby made as under:-

Thomas				<u> </u>	
Item		Description			
Code	Sutures Materials	3 2	Size	Length	Amen dments
SM- 212	Braided silk suture	3/8 Circle Cutting Needle 60 MM	1-0	76cm (±10cm)	Please read No-1 in place of 1-0 and ±15cm in place of ±10cm
SM- 216	Braided silk suture	1/2 Circle Round Body/ Taper point 25- 26 MM	3-0	45cm (± 5cm)	Please read 90 in place of 45 cm
SM- 231	Braided silk suture	Pre-cut sutures Black Braided silk sterile	1	30cm (±5cm)	Please read 75cm in place of 30 cm
SM- 255	Polyglactin/Glycoli de &lactide braided absorbable suture (dyed)	5/8 Circle Cutting 16- 17 MM	4-0	75cm (±10cm)	Please read ±15cm in place of ±10cm
SM- 269	Polyglactin/Glycolid e & lactide braided Absorbable suture Undyed	3/8 Circle Cutting 26 MM	2-0	90cm (±10cm)	Please read ±15cm in place of ±10cm
SM- 282	Polyglactin/Glycolid e &lactide Braided triclosan Coated suture	1/2 circle reverse cutting OS 40 mm	2	90cm (±10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 283	Polyglactin/Glycolid e &lactide triclosan Braided Antibacterial suture (dyed)	1/2 Circle Round Body/Taper point17- 20 MM	4-0	70cm (±10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 284	Polyglactin /Glycolide &lactide, triclosan Braided Antibacterial suture (dyed)	1/2 Circle Round Body/Taper point (Heavy) 40 MM	1	90cm (± 10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 285	Polyglactin/Glycolid e &lactide, triclosan Braided Antibacterial suture (dyed)	1/2 Circle Round Body/Taper point 20- 22 MM	3-0	70cm (± 10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 286	Polyglactin/Glycolid e &lactide, triclosan Braided Antibacterial suture (dyed)	1/2 Circle Reverse Cutting 40 MM	1	90cm (± 10cm)	Please read standard coated suture in place of triclosan Coated suture and ±15cm in place of ±10cm
SM- 294		3/8 Circle Cutting Needle 22 MM	5-0	90cm (±10cm)	Please read Cutting Needle 12-22 MM in place of





1					
	non absorbable suture				22mm and ±15cm in place of ±10cm
SM- 295	Nylon polyamide Monofilament non absorbable suture	3/8 Circle Cutting Needle	5-0	35cm (±5cm)	Please read Cutting Needle 12-22 MM
SM- 296	Nylon polyamide Monofilament non absorbable suture	3/8 Circle Cutting Needle	6-0	70cm (±10cm)	Please read Cutting Needle 12-22 MM a.nd ±15cm in place of ±10cm
SM- 299	Nylon polyamide Monofilament non absorbable suture	3/8 Circle Reverse Cutting Needle 12-13 MM	5-0	45cm (±5cm)	
8M- 316	Polyglecaprone/ Glycomer Monofilament (undyed) Suture	3/8 Circle Cutting	4-0	70cm (±10cm)	Please read Cutting Needle 12-22 MM and ±15cm in place of ±10cm
SM- 329	Polypropylene non absorbable Monofilament (dyed) suture	3/8 circle taper point 175-6, EP Double needle 8mm Polyproplyne Blue Monofilament	8-0	60cm (±10cm)	Please read 3/8 circle taper point, Double needle 8mm Polyproplyne Blue Monofilament in place of 3/8 circle taper point 175-6, EP Double needle 8mm Polyproplyne Blue Monofilament and ±15cm in place of ±10cm
8M- 330	Polypropylene non absorbable Monofilament (dyed) suture	3/8 circle taper point 175-6, EP Double needle 13mm Polyproplyne Blue Monofilament	6-0	60cm (±10cm)	3/8 circle taper point Double needle 13mm Polyproplyne Blue Monofilament in place of 3/8 circle taper point 175- 6, EP Double needle 13mm Polyproplyne Blue Monofilament and and ±15cm in place of ±10cm
SM- 331	Polypropylene non absorbable Monofilament (dyed) suture	3/8 circle taper point 175-6, EP tungsten rhenium needle BV- 175-6, Double needle 9.3 mm Polyproplyne Blue Monofilament	6-0	60c (±10cm)	3/8 circle taper point tungsten /alloy needle, Double needle 8-10.5 mm Polyproplyne Blue Monofilament in place of 3/8 circle taper point 175-6, EP tungsten rhenium needle BV- 175-6, Double needle 9.3 mm Polyproplyne Blue Monofilament. ±15cm in place of ±10cm
SM- 337	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point 25- 26 MM	3-0	90cm (±10cm)	Please read ±15cm in place of ±10cm
3M- 338	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Taper point Double Needle 26 MM	3-0	90cm (±10cm)	Please read ±15cm in place of ±10cm
3M- 339	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point Double Needle 26 MM	3-0	90cm (±10cm)	Please read ±15cm in place of ±10cm
M- 40	Polypropylene non absorbable	3/8 Circle Cutting 22- 24 MM	3-0	45cm (± 5cm)	No Change



2	T		_		
	Monofilament (dyed) suture			<u> </u>	
SM- 342	Polypropylene non absorbable Monofilament (dyed) suture	3/8 Circle Cutting 22- 24 MM	4-0	70cm (± 10cm)	Please read Cutting Needle 12-22 MM
SM- 347	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Taper Point Double Needle 26 MM	4-0	90cm (± 10cm)	Please read ±15cm in place of ±10cm
SM- 352	Polypropylene non absorbable Monofilament (dyed) suture	3/8 Circle Round Body/Taper point (282 Microns) Double Needle 11 MM	6-0	70cm (± 10cm)	Please read ±15cm in place of ±10cm and 250-350 Microns in place of 282 microns
SM- 353	Polypropylene non absorbable Monofilament (dyed) suture	3/8 Circle Round Body/Taper point, (282 Microns) Double Needle 13 MM	6-0	70cm (±10cm)	Please read ±15cm in place of ±10cm and 250-350 Microns in place of 282 microns
SM- 354	Polypropylene non absorbable Monofilament (dyed) suture	3/8 Circle Round Body/Taper point (220 Microns) Multipass BV-1 Double Needle 9-9.3 MM	6-0	70cm (±10cm)	Please read ±15cm in place of ±10cm and 200-250 Microns in place of 220 microns
SM- 358	Polypropylene non absorbable Monofilament (dyed) suture	12MM, 3/8 Circle Taper Point,	6-0	75cm (±10cm)	Please read Cutting Needle 12-22 MM Please read ±15cm in place of ±10cm
SM- 360	Polyester Braided Silicone Coated non absorbable suture	45MM, 1/2 Circle Tapercutting	5	75cm (±10cm)	Please read 45-55mm in place of 45mm Please read ±15cm in place of ±10cm
SM- 361	Polyester Braided Silicone Coated non absorbable suture	16MM, 1/2 Circle Tapercutting,	2-0	75cm (±10cm)	Please read Cutting Needle 12-22 MM Please read ±15cm in place of ±10cm
SM- 362	Polyester Braided Silicone Coated non absorbable suture	40MM, 1/2 Circle Tapercutting	2	75cm (±10cm)	Please read 40-45mm in place of 40mm Please read ±15cm in place of ±10cm
SM- 376	Knotless Sutures	26MM, 1/2 circle, taper point	3-0	15cm (±2cm)	Please read ±5cm in place of ±2cm
SM- 377	Knotless Sutures	17MM, 1/2 circle, taper point	3-0	20cm (±2cm)	Please read ±5cm in place of ±2cm
SM- 378	Knotless Sutures	27MM, 1/2 circle, taper point	2-0	15cm (± 2cm)	Please read ±5cm in place of ±2cm
SM- 383	Knotless Sutures	26MM, 1/2 circle, taper point violet triclosan coated	2-0	45cm (±5cm)	Please read standard coated in place of triclosan coated
SM- 384	Knotless Sutures	26MM, 1/2 circle, taper point violet triclosan coated	3-0	45cm (±5cm)	Please read standard coated in place of triclosan coated
SM- 385	Knotless Sutures	40MM, 1/2 circle, taper point violet triclosan coated	1-0	45cm (± 5cm)	Please read standard coated in place of triclosan coated
SM- 386	Knotless Sutures	40MM, 1/2 circle, taper point violet triclosan coated	1	45cm (± 5cm)	Please read standard coated in place of triclosan coated
SM- 391	Polytribolate coated polybutester non	13MM, 3/8 Circle Tapercutting,	5-0	75cm (±10cm)	Please read ±15cm in place of ±10cm and polybutester

A



V	.1. 1.1.				
	absorbable monofilament suture				as deleted
SM- 392	Polytribolate coated polybutester non absorbable monofilament suture	17MM, 1/2 Circle Taper Point,	4-0	90cm (±10cm)	Please read ±15cm in place of ±10cm an d polybutester as deleted
SM- 393	Polytribolate coated polybutester non absorbable monofilament suture	9MM, 3/8 Circle Taper Point	6-0	75cm (±10cm)	Please read ±15cm in place of ±10cm an d polybutester as deleted
SM- 394	Polytribolate coated polybutester non absorbable monofilament suture	8MM, 3/8 Circle Taper Point,	7-0	60cm (±10cm)	Please read ±15cm in place of ±10cm and polybutester as deleted
SM-401		Titanium ligating clips s with well defined serrati the outside clip surface, matching applier serrati interlock with clip, clip i shaped when loaded inte applier. Having lateral a transverse grooves on in surface and distal clip co providing on vessel secu	ons on ons s o the nd iner clip losure	small	Read Ligation Clip in place of Ligaclip wherever appears
SM-402	LIGACLIP	Titanium ligating clips n with well defined serration the outside clip surface, matching applier serration interlock with clip, clip is shaped when loaded into applier. Having lateral attransverse grooves on in clip surface and distal closure providing on vessecurity	ons on ons s o the nd ner	medium	Read Ligation Clip in place of Ligaclip
SM-403	LIGACLIP	Titanium ligating clips in large with well serrations on the outsi surface, matching serrations interlock wit clip is shaped when load the applier. Having late transverse grooves on clip surface and districtions on security	defined de clip applier th clip, led into ral and inner al clip vessel	medium large	Read Ligation Clip in place of Ligaclip
SM-404	LIGACLIP	Titanium ligating clips la with well defined serration the outside clip surface, matching applier serration interlock with clip, clip is shaped when loaded into applier. Having lateral and transverse ground on inner clip surface and clip closure providing on security	ons on ons othe ooves	large	Read Ligation Clip in place of Ligaclip
SM-405		Linear Cutter 55/60 mm with cartridge with 6 rows of staple line mechanism for consisten	s with		Please read 4-6 rows in place of 6 rows of staple lines



A



1		staple formation, and allowing	
		for repositioning of tissue	
		before firing. Staple row	
		extends approximately 1.5	
		staples beyond cut line. With	
		preloaded blue cartridge having	
		closed staple height of 1.5 mm.	200
SM-407	LINEAR CUTTER	Linear Cutter 75/80 mm	Please read 4-6 rows in place
		with cartridge	of 6 rows of s taple lines
		With 6 rows of staple lines with	
		mechanism for consistent	
		staple formation, and allowing	
		for repositioning of tissue	
		before firing. Staple row	
		extends approximately 1.5	
		staples beyond cut line. With	
		preloaded blue cartridge having	
		closed staple height of 1.5mm.	
SM-409	LINEAR CUTTER	Linear Cutter 100 mm with	Please read 4-6 rows in place
	DIVERNI COTTER	cartridge	of 6 rows of staple lines
		with 6 rows of staple lines with	or o rows or scapic lines
		mechanism for consistent	
		staple formation, and allowing	
	Ĭ	for repositioning of tissue	
		before firing. Staple row	
	1	extends approximately 1.5	
	T/	staples beyond cut line. With	
		preloaded blue cartridge having	
		closed staple height of 1.5 mm.	
SM-411	Endoscopic linear	(normal Shaft Length) -	Please read with articulation
	cutter	Endoscopic linear cutter	in place of single handed
		having 60 mm staple line, 57	natural articulation
		mm cut line, with single	
		handed natural articulation	
		(options of 15°, 30°, 45° on	
		either side), 360°rotation, 330-	
440		350 mm shaft length.	
SM-412	Endoscopic linear	(Bariatric Shaft Length) -	Please read with articulation
		170	
	cutter	Endoscopic linear cutter	in place of single handed
	cutter	having 60 mm staple line, 57	
	cutter	having 60 mm staple line, 57 mm cut line, with single	in place of single handed
	cutter	having 60 mm staple line, 57 mm cut line, with single handed natural articulation	in place of single handed
	cutter	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30	in place of single handed
	cutter	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side),	in place of single handed
	cutter	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm	in place of single handed
		having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length.	in place of single handed natural articulation
SM-414	Endoscopic clip	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier	in place of single handed natural articulation Please read clip length as
SM-414		having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of
SM-414	Endoscopic clip	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm.
SM-414	Endoscopic clip	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of
SM-414	Endoscopic clip	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm. Item to be added
SM-414	Endoscopic clip	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of medium large size (LT300) with	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm.
SM-414	Endoscopic clip	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of medium large size (LT300) with shaft length of 33-35mm &	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm. Item to be added
	Endoscopic clip applier	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of medium large size (LT300) with shaft length of 33-35mm & closed clip length of 8.8mm.	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm. Item to be added LT 200
SM-414 SM-415	Endoscopic clip applier Clip applier	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of medium large size (LT300) with shaft length of 33-35mm & closed clip length of 8.8mm. Multi Clip Applier	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm. Item to be added LT 200 Please read Multi clip applier
	Endoscopic clip applier	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of medium large size (LT300) with shaft length of 33-35mm & closed clip length of 8.8mm.	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm. Item to be added LT 200 Please read Multi clip applier 25-30cm with anti clip/super
	Endoscopic clip applier Clip applier	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of medium large size (LT300) with shaft length of 33-35mm & closed clip length of 8.8mm. Multi Clip Applier	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm. Item to be added LT 200 Please read Multi clip applier 25-30cm with anti clip/super interlock clip ratchet /non
	Endoscopic clip applier Clip applier	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of medium large size (LT300) with shaft length of 33-35mm & closed clip length of 8.8mm. Multi Clip Applier (Preloaded)-Open	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm. Item to be added LT 200 Please read Multi clip applier 25-30cm with anti clip/super interlock clip ratchet /non ratchet mechanism 15-30° in
	Endoscopic clip applier Clip applier	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of medium large size (LT300) with shaft length of 33-35mm & closed clip length of 8.8mm. Multi Clip Applier (Preloaded)-Open Multi clip applier 22-25cm with	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm. Item to be added LT 200 Please read Multi clip applier 25-30cm with anti clip/super interlock clip ratchet /non ratchet mechanism 15-30° in place of multi clip applier 22-
	Endoscopic clip applier Clip applier	having 60 mm staple line, 57 mm cut line, with single handed natural articulation (options of 15 degree, 30 degree, 45° on either side), 360° rotation, 430-450 mm shaft length. Multi Clip Applier (Preloaded)- Laparoscopic Endoscopic rotating clip applier 5mm with 15-20 preloaded titanium clip of medium large size (LT300) with shaft length of 33-35mm & closed clip length of 8.8mm. Multi Clip Applier (Preloaded)-Open Multi clip applier 22-25cm with anti-clip drop ratchet	in place of single handed natural articulation Please read clip length as 8.8mm to 9.1mm in place of 8.8mm. Item to be added LT 200 Please read Multi clip applier 25-30cm with anti clip/super interlock clip ratchet /non ratchet mechanism 15-30° in

W



.)		enhanced visibility, 20- 25		transverse gro ves,
		preloaded titanium clip of		30 degree
		medium size		
SM-416	Laparoscopic	Endoscopic linear cutter		Please read Eradoscopic linear
\$111- 71 0	powered linear	cartridge 35 mm having closed		cutter cartridge 30-45 mm in
	cutter 35 mm	staple height of 1.0 mm		place of Endos copic linear
	cartridge	compatible with		cutter cartridge 35 mm
	Cartriage	powered vascular stapler		attor our truge of min
		35mm.		·
SM-417	Skin stapler	Multi-directional disposable	N. Committee	Please read kryotex as deleted
		skin stapler with coated		
		rectangular staples (kryotex)		
		having dimension of 6.9mm x		
	9.	3.9mm with 35 nimbers of		
		staples with alignment		1
		indicator. Squeeze handles staple extractor.		
SM-418	Circular Stapler	Intra-luminal staplers with		Please read having controlled
		curved shaft, having controlled		tissue compression/DST
		tissue compression through		inplace of controlled tissue
		adjustable closed staple height	0	compression
		to accommodate various tissue		
		thicknesses 1.0-2.5, non-grip		
		surface, anvil detachable of		
		cartridge sizes range		
		from 21mm to 33 mm.		
SM-435	Polypropylene/	Polypropylene/ Polyester Mesh	7.6 x	Please read 7.6x15cm ±5%
SM-436	Polyester Mesh Polypropylene/	Polypropylene/ Polyester Mesh	15cm 15 x	Please read 15x15cm ±5%
9111-430	Polyester Mesh	Folypropylene/ Folyester Mesn	15 x	Flease lead 15X15cm ±5%
SM-437	Polypropylene/	Polypropylene/ Polyester Mesh	20 x	Please read 20x15cm ±5%
	Polyester Mesh	}	15cm	
SM-438	Polypropylene/	Polypropylene/ Polyester Mesh	25 x	Please read 25x20cm ±5%
	Polyester Mesh		20cm	
SM-439	Polypropylene/	Polypropylene/ Polyester Mesh	12 x	Please read12x8cm
	Polyester Mesh		8cm	±5%
SM-440	Polypropylene/	Polypropylene/ Polyester Mesh	20 x	Please read20x20cm
	Polyester Mesh		20cm	±5%
SM-441	Polypropylene/	Polypropylene/ Polyester Mesh	30 x	Please read30x20cm
	Polyester Mesh		20cm	±5%
SM-442	Tissue separator	Tissue separator mesh	12 x	Please read with/without
	mesh	{polyproplyne+ oxidized	12cm	preplaced sutures in place of
	{polyproplyne+	regenerated cellulose/		with preplaced sutures.
	oxidized	Polyester+collagen) with		Please read 12x12cm
1	regenerated	preplaced sutures		±5cm
	cellulose/ Polyester			
	collagen) with			
077 640	preplaced sutures	Wigging agreements	15	Diagram and with the
SM-443	Tissue seprator mesh	Tissue separator mesh {polyproplyne+ oxidized	15 x 10cm	Please read with/without preplaced sutures in place of
	{polyproplyne+	regenerated cellulose/	TOCH	with preplaced sutures
	oxidized	Polyester+collagen} with		Please read 15x10cm
	regenerated	preplaced sutures		±5cm
	cellulose/Polyester	F - F		
	+coll agen} with			
	preplaced			
{	Sutures		- 2	
		n A		107 Table - 0.00

A



	items in one group to co group in the tender	ompete for that	items in or tender	ne group t	o compete for that group in the	
No. 1	The Bidder must quote	at least 70% of the	4	Please the Bidder must quote at least 30% of the		
SM-45	soft polyproplyne and oxidized regenerated cellulose Inplace of	Multi layered, tissue separating mesh with large pore, soft polyproplyne and oxidized regenerated cellulose, Please Re		30.5 x 30.5cm	Please read 15x30cm ±5cm and absorbable barrier in place of oxidized regenerated cellulose	
SM-44	and oxidized regenerated cellulose	Multi layered, tissue separating mesh with large pore, soft polyproplyne and oxidized regenerated cellulose,		20cm	±5cm and absorbable barrier in place of oxidized regenerated cellulose	
SM-44	mesh{polyproplyne + oxidized regenerated cellulose/ Polyester +collagen} with preplaced sutures	Tissue separator mesh {polyproplyne+ oxidized regenerated cellulose/Polyester+collagen} with preplaced sutures		30 x 20cm	Please read with/without preplaced sutures in place of with preplaced sutures Please read 30x20cm ±5cm Please read 15x20cm	
SM-44	47 Tissue seprator mesh {polyproplyne+ oxidized regenerated cellulose/ Polyester+collagen} with preplaced sutures	Tissue seprator mesh {polyproplyne+ oxidized regenerated cellulose/Polyester+collagen} with preplaced sutures		22 x 20cm	Please read with/without preplaced sutures in place of with preplaced sutures Please read 22x20cm ±5cm	
SM-44		Tissue separator mes {polyproplyne+ oxidiz regenerated cellulose/Polyester+c with preplaced suture	ed ollagen}	20x20c m	Please read with/without preplaced sutures in place of with preplaced sutures Please read 20x20cm ±5cm	
SM-44	45 Tissue seprator mesh {polyproplyne+ oxidized regenerated cellulose/Polyester + collagen} with preplaced sutures	Tissue seprator mesh {polyproplyne+ oxidiz regenerated cellulose/Polyester+c with preplaced suture	ed ollagen}	20 x 15cm	Please read with/without preplaced sutures in place of with preplaced sutures Please read 20x15cm ±5cm	
SM-44	44 Tissue separator mesh {polyproplyne+ oxidized regenerated cellulose/Polyester +collagen} with preplaced Sutures	Tissue separator mes {polyproplyne+ oxidiz regenerated cellulose/Polyester+c with preplaced suture	ed ollagen}	15 x 15cm	Please read with/without preplaced sutures in place of with preplaced sutures Please read 15x15cm ±5cm	

Following items are added in the NIT as per the recommendations made by the technical Experts.

SM-453	Clip Applicator	Reusable Ligation Clip Applier- straight applier for small size clip 19.1cm
	Open	long
SM-454	Clip Applicator	Reusable Ligation Clip Applier- straight applier for medium size clip 19.1cm
	Open	long

SM-	Clip Applicator	Reusable Ligation Clip Applier- straight applier for Medium size clip 26.7
455	Open	cm long
SM- 456	Clip Applicator Open	Reusable Ligation Clip Applier- straight applier for large saze clip 26.7cm long
SM-	Endoscopic Clip	Endoscopic Reusable Ligation Clip Applier having alligator-ja w design, distal
457	Applicator	tip closure, rotating shaft, lateral/ transverse grooves for on vessel security, 360 rotation, for Medium size clips.
SM-	Endoscopic Clip	Endoscopic Reusable Ligation Clip Applier having alligator-ja w design, distal
458	Applicator	tip closure, rotating shaft, lateral/ transverse grooves for on vessel security, 360 rotation, for Medium large size clips with 10 mm diameter with jaw aperture 6.4mm, closed clip length 8.99 mm
SM-	Endoscopic Clip	Endoscopic Reusable Ligation Clip Applier having alligator-ja w design, distal
459	Applicator	tip closure, rotating shaft, lateral/ transverse grooves for on vessel security, 360 rotation, for large size clips.
SM-	Reusable Linear	Reusable Linear cutter/Stapler 55mm that give 57mm Staple line & 53 cut
460	Cutter Handle	line recommended for 200 clinical firings with option of u sing Knife /No Knife module to do linear cutting /stapling
SM-	Linear Cutter	Linear Cutter cartridge Blue 55 mm compatible with Reusable Linear cutter
461	Cartridge	handle 55 mm Linear Cutter cartridge Green 55 mm compatible with Reusable Linear cutter
SM- 462	Linear Cutter Cartridge	handle 55 mm
SM-	Linear Cutter	Linear Cutter knife module compatible with Reusable Linear cutter /stapler
463	Cartridge	handle 55 mm
SM-	Linear Cutter	Linear Cutter no knife module compatible with Reusable Linear cutter /stapler handle 55 mm
464 SM-	Cartridge Linear Cutter	Linear cutter with Selectable staple height 55 mm having 6 rows of staple
465	Cartridge	line, having staple height selector with cartridge having selectable closed staple height of (1.5mm/1.8mm/2mm).
SM-	Linear Cutter	Linear cutter cartridge 55 for open linear cutter compatible with selectable
466	Cartridge	staple height linear cutter 55mm. Option of closed staple height of 1.5mm/1.8mm/2mm in one cartridge only. 6 Rows 3-D staple technology.
SM-	Linear Cutter	Linear cutter with Selectable staple 75 mm having 6 rows of staple line,
467		having staple height selector with cartridge having selectable closed staple height of (1.5mm/1.8mm/2mm) & integrated new knife blade.
SM-	Linear Cutter	Linear cutter cartridge 75mm open linear cutter compatible with selectable staple height linear cutter 75mm. Option of closed staple height of
468 	Cartridge	1.5mm/1.8mm/2mm in one cartridge only. 6 Rows 3-D Staple technology.
SM-	Endoscopic	Endoscopy linear cutter 45 mm Articulating Endo Cutter with a disposable battery pack, for enhanced distal tip stability while firing, having closed
469	Linear Cutter 45mm	channel in the cartridge jaw for better stability during firing, 360, degree rotation shaft and one handed natural articulation up to 45 degrees, Precision Machined Anvil to deliver initial, system wide compression, Wide Proximal to Distal Jaw Aperture (proximal 8mm, Distal 22mm), 3 Point Gap Control for alignment and calibration throughout the 45 mm staple line, knife direction/reverse control to discontinue the firing and return the knife, interchangeable 6 row cartridge options of white, blue, gold green and black, all fits down to 12 mm trocar sleeve, 340 mm shaft length.
SM-	Endoscopic	Endoscopic linear cutter cartridge with Gripping surface Technology having 45 mm staple line, compitable with articulating endo cutter which can
470	Linear Cutter Cartridge	articulate upto 45 degree white cartridge having closed staple height of 1.0mm compatible with all range of Echelon Endoscopic linear cutter 45 mm.
SM-	Endoscopic	Endoscopic linear cutter cartridge having 45 mm staple line, compatible with
471	Linear Cutter Cartridge	articulating endo cutter which can articulate up to 45 degree blue cartridge having closed staple height of 1.5 mm compatible with all range of Echelon Endoscopic linear cutter 45 mm.
SM-	Endoscopic	Endoscopic linear cutter cartridge having 45 mm staple line, compatible with
472	Linear Cutter Cartridge	articulating endo cutter which can articulate up to 45 degree gold cartridge having closed staple height of 1.8 mm compatible with all range of Echelon Endoscopic linear cutter 45 mm.
		Landous Production of the Control of

SM-473		Endoscopic linear cutter cartridge having 45 mm staple line, compatible with
	Linear Cutter	articulating endo cutter which can articulate up to 45 degree green cartridge
	Cartridge	having closed staple height of 2.0 mm compatible with all range of Echelon Endoscopic linear cutter 45 mm.
SM-474	Endoscopic	Endoscopic linear cutter cartridge having 45 mm staple line, compatible with
	Linear Cutter	articulating endo cutter which can articulate up to 45 degree black cartridge
	Cartridge	having closed staple height of 2.3 mm compatible with all range of Echelon
		Endoscopic linear cutter 45 mm.
SM-475	1	Proximate Intra Luminal stapler with Curved shaft anvil detachable 25 mm
	Intraluminal	with powered firing, Atraumatic Gripping surface Technology with 3D
OM 476	Stapler Powered	Stapling Technology
SM-476	Intraluminal	Proximate Intra Luminal stapler with Curved shaft anvil detachable 29 mm with powered firing, Atraumatic Gripping surface Technology with 3D
	Stapler	Stapling Technology
SM-477	Powered Intra	Proximate Intra Luminal stapler with Curved shaft anvil detachable 31 mm
	luminal	should have powered firing, Atraumatic Gripping surface Technology with 3D
	Stapler	Stapling Technology
SM-478	Curved Linear	Curved Linear cutter having four curvilinear rows of staples in the width of
	Cutter	30 mm space with a 40 mm single cut line in between the staple line,
		simultaneous cutting and stapling, unique curved head design, lower pelvic access, manual pin placement option, reloadable cartridges
SM-479	Curved Linear	a) Curved blue cartridge having close staple height of 1.5mm, tactile
Din 11.5	Cutter	feedback on completion of firing sequence, new anvil, knife with every
		catridge
SM-480	Curved Linear	b) Curved, green cartridge having close staple height of 2.0 mm, tactile
	Cutter	feedback on completion of firing sequence new anvil, knife with every
SM-481	Haemorrhoidal	catridge Haemorrhoidal stapler kit consists of 33mm Haemorrhoidal Circular stapler
SW-401	Stapler	with fixed anvil
SM-482	Endocutter	60mm Articulating Endo Cutter with a disposable battery pack, for enhanced
	Powered	distal tip stability while firing, having closed channel in the cartridge jaw for
		better stability during firing, 360 degree rotation shaft and one handed
		natural articulation up to 45 degrees, Precision Machined Anvil to deliver
	!	initial, system wide compression, Wide Proximal to Distal Jaw Aperture (proximal 8mm, Distal 22mm), 3 Point Gap Control for alignment and
		calibration throughout the 60mm staple line, knife direction/reverse control
		to discontinue the firing and return the Knife, interchangeable 6 row
		cartridge options of closed staple height of 1.0, 1.5,1.8, 2.0 and 2.3 mm, all
		fits down to 12 mm trocar sleeve, 440 mm shaft length.
SM-483	Endocutter	Endoscopic linear cutter cartridge with Gripping Surface Technology having
	Powered Cartridge	60mm staple line, compatible with articulating upto 45° white cartridge having closed staple height of 1.0,1.5mm, 1.8mm,2.0mm, 2.3mm compatible
ĺ	Cartridge	with all range of Echelon Endoscopic linear cutter 60mm.
SM-484	Intracorpeal	Fibrin Sealant alongwith Dual applicator for open and Laparoscopic surgery
	Tissue Glues	both Rigid (30-35cm) and flexible (40-45cm)
	for Mesh	G .
	Fixation	
SM-485		antibacterial barb suture with fixation tab 1 No Rev Cut, 36MM, 45cm
SM-486		antibacterial barb suture with fixation tab 1 No Rev Cut, 48MM, 45cm
SM-487		antibacterial barb suture with fixation tab 2/0 Rev Cut, 26MM, 45cm
SM-488		antibacterial barb suture with LOOP FLEXIBLE 3/0 RC, 26MM, 45cm
SM-489 SM-490		antibacterial barb suture with LOOP FLEXIBLE 1/0 30CM 36MM
SM-490	36MM, taperpo	intibacterial barb suture with fixation TAB 1/0, 45CM, 1/2CIRC,
SM-491		ntibacterial barb suture with fixation TAB 1No, 1/2circ, 36mm, 45cm
SM-492		ntibacterial 1 No ½ CIRC, Taperpoint, 40mm, 90cm
SM-493	 	ntibacterial 1/0 No ½ CIRC, Taperpoint, 48mm, 90cm
SM-494		ntibacterial 2/0 No ½ CIRC, Taperpoint, 26mm, 70cm
SM-495		ntibacterial 3/0 No ½ CIRC, Taperpoint, 17mm, 70cm
		ntibacterial 4/0 No ½ CIRC, Taperpoint, 17mm, 70cm





SM-497	Polydixanone 5/0 3/8 CIRC, 17MM DOUBLE NEEDLE, 70cm
SM-498	Polydixanone6/0 3/8 CIRC, Taperpoint, 11MM 45cm
SM-499	Polydixanone antibacterial 6/0 3/8 CIRC, Taperpoint, 8MM, 45cm
SM-500	Polydixanone Port Closure NEEDLE 22mm cutting
SM-501	FAST Absorbable Suture (Glycolide 90% with 10% L-Lactide as two components with
	Standard Surface Coating), colour violet 1/0 TAPERPOINT 40MM, 110C M, TAPER
	POINT ½ CIRC,
SM-502	FAST Absorbable Suture (Glycolide 90% with 10% L-Lactide as two components with
	Standard Surface Coating), colour violet 1/0 TAPERPOINT 31MM, 110C M, 1/2CIRC
SM-503	FAST Absorbable Suture (Glycolide 90% with 10% L-Lactide as two components with
	Standard Surface Coating), colour violet 3/0 TAPERPOINT 26MM, ½ CIRC, 70CM
SM-504	FAST Absorbable Suture (Glycolide 90% with 10% L-Lactide as two components with
	Standard Surface Coating), colour violet 2/0 TAPER CUT 31MM, 110CM, 1/2CIRC
SM-505	PLOYGLYPRONE 26MM CUTTING FOR SUBCUTE, 3/0, 3/8 circ reverse cutting 70cm
SM-506	BIOPATCH Protective Disc with CHG – BIOPATCH ANTIMICROBIAL DRESSING, 1
	disk (2.5cm) w/4.0mm centre hole
SM-507	BIOPATCH Protective Disc with CHG – BIOPATCH ANTIMICROBIAL DRESSING, 3/4
	disk (1.9cm) w/1.5mm centre hole
SM-508	BIOPATCH protective Disc with CHG - BIOPATCH ANTIMICROBIAL DRESSING, 1
	disk (2.5cm) w/7.0mm centre hole
SM-509	Oxidized Regenerated Cellulose adhesion barrier – 3x4" inch
SM-510	Partially Absorbable Mesh Made of Polyglycaprone & Polypropylene 15X15
SM-511	Partially Absorbable Mesh Made of Polyglycaprone & Polypropylene 30X3
SM-512	Partially Absorbable Mesh Made of Polyglycaprone & Polypropylene 10X15

In addition to above please read:

1. Only technical/financial bids uploaded as per revised corrigendum shall be considered for evaluation. All the bidders are requested to keep themselves updated and submit their e- bids through e-portal www.jktenders.gov.in as per specifications & BOQs. The amendments / modifications shall be available on www.jktenders.gov.in, and www.jkmsclbusiness.com.

No. JKMSCL/Suture/Corg/450/2021/ 17045-47

Copy for information to the:-

1. Managing Director, JKMSCL

General Manager (Adm), JKMSCL.
 Assistant programmer to uploaded on www.jktenders.gov.in and on www.jkmsclbusiness.com

4. Office Record.