



# JAMMU AND KASHMIR MEDICAL SUPPLIES CORPORATION LTD.

(Public Sector Undertaking of the Government of Jammu and Kashmir)

Corporate Head Office: Plot No. 58, Friends Colony Satyam Road Trikuta Nagar Jammu

Corporate Office: Opp J&K Motor Garage Deptt near Hajj House Bemina Srinagar

Telefax: 0191-2478842; 0194-2490662 email: gmj@jkmscl@gmail.com; website: www.jkmsclbusiness.com

## C O R R I G E N D U M

In continuation to e-NIT No: NIT No: JKMSCL/SUTURE MATERIAL/686 Dated: 01-04-2026 invited for procurement of **Suture material**, following amendments are hereby made:

1. Sub-categorization of sutures into Vital sutures (Group A1), Minor Sutures (Group A2) and miscellaneous (Group A3) are uploaded as separate tender items, replacing the existing list of sutures.
2. Similarly, Stapling Devices, Hemostats and Meshes are categorized as Group B, C and D.

Details are as under:

### Group A1 Vital Sutures

S No.	Item Code	Name of the Item	Technical Specification		
			Needle Type & Size	Size	Length
	<b>Polyglactin/ Glycolide &amp; lactide braided absorbable sutures (dyed)</b>	<b>Sutures Materials</b>			
1	SM- 253	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/2 Circle Round Body/Taper point 20 MM	4-0	70cm ( ± 10cm)
2	SM-254	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	5/8 Circle Cutting 16-17 MM	4-0	45cm ( ± 5cm)
3	SM-255	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	5/8 Circle Cutting 16-17 MM	4-0	75 cm ( ± 15cm)
4	SM- 256	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	3/8 Circle Reverse Cutting needle 12 MM	5-0	45cm ( ± 5cm)
5	SM-257	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	3/8 Circle Cutting needle 12 MM	5-0	45cm ( ± 5cm)
6	SM- 258	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/4 Circle Spatulated MICRO-POINT Double Needle 8 MM	6-0	30cm ( ± 5cm)

W

LA

7	SM-259	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/4 Circle Spatulated MICRO-POINT Double Needle 8 MM	6-0	45cm (± 5cm)
8	SM-260	Polyglactin/Glycolide & lactide braided Absorbable suture (Rapid absorption) Undyed	3/8 Circle Reverse Cutting 24-26 MM	3-0	75cm (± 10cm)
9	SM-261	Polyglactin/Glycolide & lactide braided Absorbable suture (Rapid absorption) Undyed	3/8 Circle Cutting , needle 12-16 MM	4-0	45cm (± 5cm)
10	SM-262	Polyglactin braided Absorbable suture (Rapid absorption) Undyed	3/8 Circle Cutting , needle 16 MM	4-0	75cm (± 10cm)
11	SM-263	Polyglactin/Glycolide & lactide braided Absorbable suture (Rapid absorption) Undyed	3/8 Circle Cutting , needle 16-19 MM	3-0	45cm (± 5cm)
12	SM-264	Polyglactin/Glycolide & lactide braided Absorbable suture (Rapid absorption) Undyed	3/8 Circle Cutting , needle 16-19 MM	3-0	75cm (± 10cm)
13	SM-265	Polyglactin/Glycolide & lactide braided Absorbable suture (Rapid absorption) Undyed	1/2 Circle Cutting, 22-26 MM	3-0	45cm (± 5cm)
14	SM-266	Polyglactin/Glycolide & lactide braided Absorbable suture (Rapid absorption) Undyed	1/2 Circle Cutting , 22-26 MM	3-0	75cm (± 10cm)
15	SM-267	Polyglactin/Glycolide & lactide braided Absorbable suture (Rapid absorption) Undyed	3/8 Circle Cutting , needle 12-16 MM	4-0	70cm (± 10cm)
16	SM-268	Polyglactin/Glycolide & lactide braided Absorbable suture (Rapid absorption) Undyed	3/8 Circle Reverse Cutting 11 MM	5-0	45cm (± 5cm)
17	SM-269	Polyglactin/Glycolide & lactide braided Absorbable suture Undyed	3/8 Circle Cutting 26 MM	2-0	90cm (± 15cm)
18	SM-270	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/2 Circle Tapercut , 36 MM	0	90cm (± 10cm)
19	SM-271	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/2 Circle Round Body/Taper point 30 MM	0	80cm (± 10cm)
20	SM-272	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/2 Circle Reverse Cutting OS 40 MM	1	90cm (± 10cm)

W

SA

21	SM- 273	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/2 Circle Round Body/Taper point (Heavy) 40 MM	1	90cm ( ± 10cm)
22	SM- 274	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/2 Circle Round Body/Taper point 40 MM	1-0	70cm ( ± 10cm)
23	SM- 275	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/2 Circle Round Body/Taper point 30 MM	2-0	90cm ( ± 10cm)
24	SM-276	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	3/8 Circle Cutting 26 MM	2-0	90cm ( ± 10cm)
25	SM-277	Polyglactin/Glycolide & lactide braided suture	3/8 Circle Cutting 22-24 MM	3-0	45cm ( ± 5cm)
26	SM-278	Polyglactin/Glycolide & lactide braided suture (dyed) Absorbable	3/8 Circle Cutting 22-24 MM (±10%) 3-0 45cm ( ± 5cm)	3-0	45cm ( ± 10cm)
27	SM- 279	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/2 Circle Round Body/Taper point 30 MM	2-0	70cm ( ± 10cm)
28	SM- 280	Polyglactin/Glycolide & lactide braided suture (dyed) Absorbable	1/2 Circle Round Body/Taper point 20-22 MM	3-0	70cm ( ± 10cm)
29	SM-281	Polyglactin/Glycolide & lactide braided absorbable suture (dyed)	1/2 Circle Cutting 22-24 MM	3-0	80cm ( ± 10cm)
30	SM- 282	Polyglactin/Glycolide & lactide Braided Standard Quoted Suture	1/2 circle reverse cutting OS 40 mm	2	90cm ( ± 15cm)
31	SM- 283	Polyglactin/Glycolide & lactide Standard Quoted Suture Braided suture (dyed)	1/2 Circle Round Body/Taper point 17-20 MM	4-0	70cm ( ± 15cm)
32	SM-284	Polyglactin/Glycolide & lactide Standard Quoted Suture Braided standard quoted suture (dyed)	1/2 Circle Round Body/Taper point (Heavy) 40 MM	1	90cm ( ± 15cm)
33	SM-285	Polyglactin/Glycolide & lactide Standard Quoted Suture Braided standard quoted suture (dyed)	1/2 Circle Round Body/Taper point 20-22 MM	3-0	70cm ( ± 15cm)
34	SM-286	Polyglactin/Glycolide & lactide Standard Quoted Suture Braided standard quoted suture (dyed)	1/2 Circle Reverse Cutting 40 MM	1	90cm ( ± 15cm)
35	SM- 287	Polyglactin/Glycolide & lactide braided suture (dyed) Absorbable	3/8 Circle Spatulated MICRO-POINT needle, 30 cm	7-0	30cm ( ± 5cm)

W

ⓧ

36	SM- 288	Polyglactin/Glycolide &lactide braided suture (dyed) Absorbable	Round Body/Taper point needle, 3/8 circle Needle length:6.50 mm	8-0	30 cm ( ± 5cm)
37	SM-289	Polyglactin/Glycolide &lactide braided suture (dyed) Absorbable	cutting needle,3/8 circle spatulated Needle length 6 mm	8-0	30 cm ( ± 5cm)
38	SM- 504	Polyglactin/Glycolide &lactide braided suture (dyed) Absorbable	Fast Absorbable Suture ( Glycolide 90% with 10% L- lactide as two components with standard surface coating), Colour Violet, 2/0 Tapercut 31-37 mm, 1/2 CIRC	2-0	100cm ( ± 10cm)

S No.	Item Code	Name of the Item	Technical Specification		
			Needle Type & Size	Size	Length
	<b>Polyglecaprone / Glycomer Monofilament (undyed) Sutures</b>	<b>Sutures Materials</b>			
39	SM- 314	Polyglecaprone/ Glycomer Monofilament (undyed) Suture	1/2 Circle Oval Round Body/Taper point, 22 MM	3-0	70cm (±10cm)
40	SM- 315	Polyglecaprone/ Glycomer Monofilament (undyed) Suture	3/8 Circle Cutting 25 MM	3-0	70cm (±10cm)
41	SM- 316	Polyglecaprone/ Glycomer Monofilament (undyed) Suture	3/8 Circle reverse Cutting needle 12-22 MM	4-0	70cm (±15cm)
42	SM- 317	Polyglecaprone/ Glycomer Monofilament (undyed) Suture	3/8 Circle Cutting 13-16 MM	5-0	45cm (± 5cm)

S No.	Item Code	Name of the Item	Technical Specification		
			Needle Type & Size	Size	Length
	<b>Polydioxanone/Polyglyconate monofilament sutures</b>	<b>Sutures Materials</b>			
43	SM-320	polydioxanone/Polyglyconate monofilament suture	1/2 Circle Round Body/Taper point CTX 48 MM	1	150cm ( ± 10cm)
44	SM- 321	polydioxanone/Polyglyconate monofilament Suture	1/2 Circle Round Body/Taper point Double Needle 13 MM	5-0	70cm ( ± 10cm)

WV.

DA

45	SM- 322	polydioxanone/Polyglyconate monofilament Suture	1/2 Circle, Round Body/Taper point, double needle, 17 mm	4-0	90cm ( ± 10cm)
46	SM- 324	polydioxanone/Polyglyconate monofilament suture	1/2 Circle, Round Body/Taper point, 20mm	3-0	70cm ( ± 10cm)
47	SM- 325	polydioxanone/Polyglyconate monofilament suture	Round bodied, 1/2 circle, 13mm	4-0	45cm ( ± 5cm)
48	SM- 326	polydioxanone/Polyglyconate monofilament suture	Round bodied, 1/2 circle, 20mm	4-0	70cm ( ± 10cm)
49	SM- 327	polydioxanone/Polyglyconate monofilament suture 1/2 Circle Round Body/Taper Point 13 MM(±10%) 5-0 70cm ( ± 10cm)	polydioxanone/Polyglyconate monofilament suture Round bodied, 13mm 1/2 Circle 70cm ( ± 10cm)	5-0	70cm (±10cm)
50	SM- 492	POLYDIXANONE/POLYGLYCONATE 1 no. 1/2 circle taperpoint 40mm 90cm	POLYDIXANONE/POLYGLYCONATE 1 no. 1/2 circle taperpoint 40mm 90cm	1	90cm ( ± 10cm)
51	SM- 494	POLYDIXANONE/POLYGLYCONATE 2/0 . 1/2 circle taper-point 26mm 90cm	POLYDIXANONE/POLYGLYCONATE 2/0 . 1/2 circle taper-point 26mm 90cm	2/0	70cm (±10cm)
52	SM- 495	POLYDIXANONE/POLYGLYCONATE 3/0 . 1/2 circle taper-point 26mm 90cm	POLYDIXANONE/POLYGLYCONATE 3/0 . 1/2 circle taper-point 26mm 90cm	3/0	70cm (±10cm)
53	SM-496	polydioxanone/Polyglyconate monofilament suture	POLYDIXANONE/POLYGLYCONATE 4/0 . 1/2 circle taperpoint 17mm 75cm	4/0	70cm (±10cm)
54	SM- 497	POLYDIXANONE/POLYGLYCONATE 5/0 . 3/8 circle taperpoint 17mm Double arm 70cm	POLYDIXANONE/POLYGLYCONATE 5/0 . 3/8 circle taperpoint 17mm Double arm 70cm	5/0	70cm (±10cm)
55	SM- 498	POLYDIXANONE/POLYGLYCONATE 6/0 . 3/8 circle taperpoint 11mm Double arm 45-75cm	POLYDIXANONE/POLYGLYCONATE 6/0 . 3/8 circle taperpoint 11-13mm Double arm 45-75cm	6/0	45cm ( ± 5cm)
56	SM- 499	POLYDIXANONE/POLYGLYCONATE 6/0 . 3/8 circle taperpoint 8mm Double arm 45-75cm	POLYDIXANONE/POLYGLYCONATE 6/0 . 3/8 circle taperpoint 8mm Double arm 45-75cm	6/0	45cm ( ± 5cm)

WV

S No.	Item Code	Name of the Item	Technical Specification		
			Needle Type & Size	Size	Length
	<b>Polypropylene non absorbable Monofilament (dyed sutures</b>	<b>Sutures Materials</b>			
57	SM- 329	Polypropylene non absorbable Monofilament (dyed) suture	3/8 circle taper point , Double needle 8mm Polypropylene Blue Monofilament	8-0	60cm (± 15cm)
58	SM- 330	Polypropylene non absorbable Monofilament (dyed) suture	3/8 circle taper point Double needle 13mm Polypropylene Blue Monofilament	6-0	60cm (± 15cm)
59	SM- 331	Polypropylene non absorbable Monofilament (dyed) suture	3/8 circle taper point tungsten /Alloy needle , Double needle, 8-10.5 mm Polypropylene Blue Monofilament	6-0	60cm (± 15cm)
60	SM- 332	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point 27-30 MM	0	80cm (± 10cm)
61	SM- 334	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point (Heavy) 27-30 MM	1	80cm (± 10cm)
62	SM- 335	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point 26-30 MM	2-0	85cm (± 10cm)
63	SM- 336	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Tapercut V- 5 Double Needle 17 MM	2-0	90cm (± 10cm)
64	SM- 337	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point 25-26 MM	3-0	90cm (± 15cm)
65	SM- 338	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Taper point Double Needle 26 MM	3-0	90cm (± 15cm)
66	SM-339	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point Double Needle 26 MM	3-0	90cm (± 15cm)
67	SM- 340	Polypropylene non absorbable	3/8 Circle Cutting 22- 24 MM	3-0	45cm (± 5cm)

✓

✓

		Monofilament (dyed) suture			
68	SM-341	Polypropylene non absorbable Monofilament (dyed) suture	3/8 Circle Cutting 22- 24 MM	3-0	70cm ( ± 10cm)
69	SM- 342	Polypropylene non absorbable Monofilament (dyed) suture	3/8 Circle Cutting 12- 22 MM	4-0	70cm ( ± 10cm)
70	SM- 343	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point Double Needle 17 MM	4-0	90cm ( ± 10cm)
71	SM- 344	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point 16- 17mm	4-0	80cm ( ± 10cm)
72	SM-348	Polypropylene non absorbable Monofilament (dyed) suture	1/2 Circle Round Body/Taper point, Double Needle 13 MM	5-0	75cm ( ± 10cm)
73	SM-349	Polypropylene non absorbable Monofilament (dyed) suture	3/8 Circle Round Body/Taper point Double Needle 12 MM	5-0	70cm ( ± 10cm)
74	SM- 351	Polypropylene non absorbable Monofilament (dyed) suture	3/8 Circle Round Body/Taper point Micro Point Double Needle 8 MM	6-0	70cm ( ± 10cm)
75	SM-352	Polypropylene non absorbable Monofilament (dyed) suture	Polypropylene non absorbable Monofilament (dyed) suture 3/8 Circle Round Body/Taper point (250-350 Microns) Double Needle 11 MM 80cm ( ± 15cm)	6-0	70cm ( ± 15cm)
76	SM- 353	Polypropylene non absorbable Monofilament (dyed) suture	3/8 Circle Round Body/Taper point, (250-350 Microns) Double Needle 13 MM	6-0	70cm ( ± 15cm)
77	SM- 354	Polypropylene non absorbable Monofilament (dyed) suture 3/8 Circle Round Body/Taper Point Double Needle 9- 9.3 MM(±10%) 6-0 70cm ( ± 10cm)	Polypropylene non absorbable Monofilament (dyed) suture 3/8 Circle Round Body Double Needle 9-9.3 mm ( ± 10%) 6-0 70cm ( ± 10cm)	6-0	70cm ( ± 15cm)

W

J

78	SM- 355	Polypropylene non absorbable Monofilament (dyed) Suture	3/8 Circle Taper Point , Double Needle 8 MM	7-0	60cm ( ± 10cm)
79	SM-356	Polypropylene non absorbable Monofilament (dyed) Suture	3/8 circle, taper point, Everpoint, Double needle 8 mm	7-0	80cm ( ± 10cm)

S No.	Item Code	Name of the Item	Technical Specification		
			Needle Type & Size	Size	Length
		<b>Sutures Materials</b>			
		<b>Polyester Braided Silicone Coated non absorbable sutures</b>			
80	SM- 360	Polyester Braided Silicone Coated non absorbable suture	45-55 mm, 1/2 Circle Taper cutting	5	75cm ( ± 15cm)
81	SM- 361	Polyester Braided Silicone Coated non absorbable suture	12-22 mm, 1/2 Circle Tapercut, Green braided	2-0	75cm ( ± 15cm)
82	SM- 362	Polyester Braided Silicone Coated non absorbable suture	40-45mm, 1/2 Circle Tapercut	2	75cm ( ± 15cm)
83	SM- 363	Polyester Braided Silicone Coated non absorbable suture	1/2 Circle Tapercut V- 5 Double Needle 16-17 MM Polybutylate coated braided polyester Excel WHITE Braided 90 cm Firm PTFE POLYMER PLEDGETS 7mm x 3mm x 1.5mm	2-0	90cm ( ± 10cm)
84	SM- 364	Polyester Braided Silicone Coated non absorbable suture	1/2 Circle Tapercut V- 5 Double Needle 17 MM Polybutylate coated braided polyester Excel Green Braided 90 cm Firm PTFE POLYMER PLEDGETS 7mm x 3mm x 1.5mm (10 Suture - Per Foil)	2-0	90cm ( ± 10cm)
85	SM- 365	Polyester Braided Silicone Coated non absorbable suture	1/2 Circle Tapercut V- 7 Double Needle 26 MM Polybutylate coated braided polyester Excel Green Braided 90 cm Firm PTFE POLYMER	2-0	90cm( ± 10cm)

W

			PLEDGETS 7mm x 3mm x 1.5mm		
86	SM- 366	Polyester Braided Silicone Coated non absorbable suture	1/2 Circle Tapercut V- 7 Double Needle 26 MM Polybutylate coated braided polyester Excel Green Braided 90 cm Firm PTFE POLYMER PLEDGETS 7mm x 3mm x 1.5mm	2-0	90cm ( ± 10cm)
87	SM- 368	Polyester Braided Silicone Coated non absorbable suture	1/2 Circle Tapercut 2 x V-5 Double Needle 17 MM Polybutylate coated braided polyester Excel Green Braided	2-0	90cm ( ± 10cm)
88	SM- 370	Polyester Braided Silicone Coated non absorbable suture	1/2 Circle Tapercut Double Needle 25-26 MMPolybutylate coated braided polyester Excel White Braided	2-0	90cm ( ± 10cm)
89	SM- 371	Polyester Braided Silicone Coated non absorbable suture	1/2 Circle Round Body/Taper point SH Needle Needle 26 MM Polybutylate coated braided polyester EXCEL Green Braided	3-0	75cm ( ± 10cm)
90	SM- 372	Polyester Braided Silicone Coated non absorbable suture	1/2 Circle Round Body/Taper point SH Double Needle 26-35 MM Polybutylate coated braided polyester Excel Green Braided	3-0	90cm ( ± 10cm)
91	SM-375	Polyester Braided Silicone Coated non absorbable suture	1/2 circle needle 26 mm, Round Body/Taper point, green braided	2-0	75cm ( ± 10cm)

S No.	Item Code	Name of the Item	Technical Specification		
			Sutures Materials	Needle Type & Size	Size
92	SM- 376	Knotless Sutures	26MM, 1/2 circle, taper point.	3-0	15cm ( ± 5cm)
93	SM- 377	Knotless Sutures	17MM, 1/2 circle, taper point	3-0	15cm ( ± 5cm)
94	SM- 378	Knotless Sutures	36MM, 1/2 circle, taper point	2-0	15cm ( ± 2cm)

W

✓

95	SM-379	Knotless Sutures	17MM, 1/2 circle, taper point	2-0	20cm ( ± 2cm)
96	SM- 381	Knotless Sutures	36mm, 1/2 Circle reverse cutting	1	24cm ( ± 2cm)
97	SM- 383	Knotless Sutures	26MM 1/2 TAPER POINT WITH STANDARD COATING 2-0 45CM	2-0	45cm ( ± 5cm)
98	SM- 384	Knotless Sutures	26MM 1/2 TAPER POINT WITH STANDARD COATING 3-0 45CM	3-0	45cm ( ± 5cm)
99	SM- 385	Knotless Sutures	40MM 1/2 TAPER POINT WITH STANDARD COATING 1-0 45CM	1-0	45cm ( ± 5cm)
100	SM-386	Knotless Sutures	40MM 1/2 Circle taper point with Standard Coated 1 NO. 45CM (± 5 CM)	1	45cm ( ± 5cm)
101	SM- 485	Knotless Sutures	KNOTLESS SUTURE POLYDIXANONE/POLYGLYCONATE symmetric with fixation tab/loop size 1NO.Reverse cutting 36mm, 45CM (±5 CM)	1	45cm ( ± 5cm)
102	SM-486	Knotless Sutures	POLYDIXANONE/ POLYGLYCONATE symmetric with fixation tab/loop size 1 NO. Taper point 48mm, 45CM (±5 CM)	1	45cm ( ± 5cm)
103	SM- 487	Knotless Sutures	KNOTLESS SUTURE POLYDIXANONE/POLYGLYCONATE symmetric with fixation tab/loop size 2/0 taper point 26mm, 45CM (±5 CM)	2/0	45cm ( ± 5cm)
104	SM- 488	Knotless Sutures	Polydixanone/polyglyconate Barb suture with loop Flexible 3/0 Taper Point, 26mm, 45cm	3/0	45cm ( ± 5cm)
105	SM- 489	Knotless Sutures	Polydixanone Barb suture with loop Flexible 1/0 30 cm, 36mm	1/0	30cm (±2cm)
106	SM- 490	Knotless Sutures	Polydixanone Barb suture symmetric with fixation Tab 1/0, 45 cm, 1/2CIRC, 36 mm, taperpoint	1/0	45cm (±5cm)
107	SM-776	Knotless Barb Sutures	Unidirectional Barbed, monofilament, Polybutester non absorbable suture ½ circle, Taper Point/ Round Body Needle 37mm	1	45cm (±5cm)

✓

✓

108	SM-777	Knotless Barb Sutures	Unidirectional Barbed, monofilament, Polybutester non absorbable suture ½ circle, Taper Point/ Round Body Needle 37mm	1-0	45cm (±5cm)
109	SM-778	Knotless Barb Sutures	Unidirectional Barbed, monofilament, Polybutester non absorbable suture ½ circle, Taper Point/ Round Body Needle 37mm	1-0	15cm (±5cm)

<b>Group A2 Trivia/Minor Sutures</b>					
<b>Sr No</b>	<b>Item Code</b>	<b>Sutures Materials</b>	<b>Needle Type &amp; Size</b>	<b>Size</b>	<b>Length</b>
1	SM-209	Braided silk suture	3/8 Circle Reverse Cutting 45 MM	0	76cm (± 10cm)
2	SM-211	Braided Silk Suture	1/2 Circle Round Body/ Taper point 30 MM	0	80cm (±10cm)
3	SM-212	Braided Silk Suture	3/8 Circle Cutting Needle 60 MM	1	76cm (± 5cm)
4	SM-213	Braided Silk Suture	1/2 Circle Round Body/ Taper point 30 MM	2-0	80cm (±10cm)
5	SM-214	Braided silk suture	3/8 Circle Reverse Cutting 39- 45 MM	Feb-00	75cm (± 10cm)
6	SM-215	Braided Silk Suture	1/2 Circle Round Body/ Taper point 20 MM	3-0	75cm (±10cm)
7	SM-216	Braided silk suture	Braided Silk Suture 1/2 Circle Round Body/ Taper point 25- 26 MM 90cm ( ± 15cm)	3-0	90cm ( ± 15cm)
8	SM-218	Braided silk suture	3/8 Circle Reverse Cutting 26 MM	Mar-00	76cm ( ± 10cm)
9	SM-219	Braided silk suture	3/8 Circle Cutting Needle 16 MM	Mar-00	90cm ( ± 10cm)
10	SM-220	Braided silk suture	3/8 Circle Cutting Needle 16 MM	4-0	90cm ( ± 10cm)
11	SM-221	Braided silk suture	3/8 Circle Reverse Cutting Needle, 12 MM-16mm	5-0	90cm ( ± 10cm)
12	SM-222	Braided silk suture	3/8 Circle Round Body/Taper point 10-12 MM 76cm ( ± 10cm)	5-0	76cm ( ± 10cm)
13	SM-223	Braided silk suture	Black Braided silk	4-0	25mtr ( ± 3mtr)
14	SM-224	Braided silk suture	Black Braided silk	3-0	25mtr ( ± 3mtr)
15	SM-225	Braided silk suture	Black Braided silk	2-0	25mtr ( ± 3mtr)
16	SM-226	Braided silk suture	Black Braided silk	0	25mtr ( ± 3mtr)

SV

17	SM-227	Braided silk suture	Black Braided silk	1	25mtr ( ± 3mtr)
18	SM-228	Braided silk suture	Pre-cut sutures Black Braided silk sterile	0	35cm ( ± 5cm)
19	SM-230	Braided silk suture	Pre-cut sutures Black Braided silk sterile	1	75cm ( ± 10cm)
20	SM-234	Braided silk suture	Pre-cut sutures Black Braided silk sterile	2-0	75cm ( ± 10cm)
21	SM-235	Braided silk suture	Pre-cut sutures Black Braided silk sterile	2-0	35cm ( ± 5cm)
22	SM-236	Braided silk suture	Pre-cut sutures Black Braided silk sterile	3-0	35cm ( ± 5cm)
23	SM-238	Braided silk suture	Pre-cut sutures Black Braided silk sterile 60-75cm	4-0	60-75cm
24	SM-240	Braided silk suture	3/8 Circle micropoint Spatulated 6 MM Virgin Silk	8-0	38cm ( ± 5cm)

S No.	Item Code	Name of the Item	Technical Specification		
			Needle Type & Size	Size	Length
	<b>Surgical Gut/ Synthetic Polyglytone Sutures</b>	<b>Sutures Materials</b>			
25	SM-245	Surgical gut/ Synthetic Polyglytone Suture	3/8 Circle Round Body/Tapper point 12-13 MM (±10%) Size 1 No.	5-0	76cm ( ± 10cm)
26	SM-246	Surgical gut/Synthetic Polyglytone Suture	1/2 Circle Round Body/ Taper point (Heavy) 40 MM	1	76cm ( ± 10cm)
27	SM-249	Surgical gut/Synthetic Polyglytone Suture	1/2 Circle Round Body/Taper point 26 MM	2-0	76cm ( ± 10cm)
28	SM-251	Surgical gut/ Synthetic Polyglytone Suture	1/2 Circle Round Body/Taper point 20 MM	3-0	76cm ( ± 10cm)

S No.	Item Code	Name of the Item	Technical Specification		
			Needle Type & Size	Size	Length
	<b>Nylon polyamide Monofilament non absorbable sutures</b>	<b>Sutures Materials</b>			
29	SM-290	Nylon polyamide Monofilament non absorbable Suture	Nylon polyamide Monofilament non absorbable Suture	0	150cm (±10cm)

W

J

			loop 1/2 Circle Round Body/Taper point 40 MM 150cm (±10cm) 1/0		
30	SM- 291	Nylon polyamide Monofilament non absorbable suture 1/2 Circle Round Body (Heavy) 48 MM(±10%) 1 150cm ( ± 10cm)	Nylon polyamide Monofilament non absorbable suture loop 1/2 Circle Round Body/Tapper Point (Heavy) 48 MM-50mm 150cm (±10cm)	1	150cm (±10cm)
31	SM- 292	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Reverse Cutting 39-45 MM	2-0	70cm ( ± 10cm)
32	SM-293	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Cutting Needle 16-19 MM	3-0	80cm ( ± 10cm)
33	SM- 294	Nylon polyamide Monofilament non absorbable suture	3/8 Circle Cutting Needle 12-22 MM, Hemoseal	5-0	90cm ( ± 15cm)
34	SM-295	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Cutting Needle size 12-22mm	5-0	40cm ( ± 5cm)
35	SM- 296	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Cutting Needle 12-22mm	6-0	70cm ( ± 15cm)
36	SM- 297	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Reverse Cutting 26 MM	3-0	70cm ( ± 10cm)
37	SM- 298	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Cutting Needle 16-19 MM	4-0	90cm ( ± 10cm)
38	SM- 299	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Reverse Cutting Needle 12-13 MM	5-0	45cm ( ± 10cm)
39	SM-300	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Reverse Cutting Needle 12-13 MM	5-0	90cm ( ± 10cm)
40	SM- 301	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Round Body/Taper point Micro Point Double Needle 4 MM	10-0	5cm
41	SM-302	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Round Body/Taper point Micro Point 4 MM	10-0	25cm ( ± 5cm)
42	SM- 303	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Round Body/Taper point Micro Point 6 MM	8-0	38cm ( ± 5cm)
43	SM- 304	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Spatulated Micro Point Double Needle 6 MM	10-0	38cm ( ± 5cm)

W

SA

44	SM-305	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Reverse Cutting 39-45 MM	2-0	70cm ( ± 10cm)
45	SM-306	Nylon polyamide Monofilament non absorbable Suture	3/8 Circle Reverse Cutting 24-26 MM	3-0	70cm ( ± 10cm)
46	SM-307	Nylon polyamide Monofilament non absorbable Suture	Straight Cutting Needle 60 MM	0	75cm ( ± 10cm)
47	SM-308	Nylon polyamide Monofilament non absorbable Suture	Nylon pre cut ties,	1	75cm ( ± 10cm)
48	SM-309	Nylon polyamide Monofilament non absorbable Suture	Nylon pre cut ties	2-0	75cm ( ± 10cm)
49	SM-310	Nylon polyamide Monofilament non absorbable Suture	Nylon pre cut ties	3-0	75cm ( ± 10cm)
50	SM-311	Nylon polyamide Monofilament non absorbable Suture	Nylon pre cut ties	4-0	45cm ( ± 5cm)
51	SM- 312	Nylon polyamide Monofilament non absorbable Suture	3/8 circle. Needle length: 11.00 mm.	7-0	45cm ( ± 5cm)
52	SM- 313	Nylon polyamide Monofilament non absorbable Suture	3/8 circle. Needle length: 6.50 mm.	9-0	30cm ( ± 5cm)

<b>Group A3 Miscellaneous</b>					
<b>S No</b>	<b>Item Code</b>	<b>Name of the Item</b>	<b>Technical Specification</b>		
	<b>Miscellaneous</b>	<b>Sutures Materials</b>	<b>Needle Type &amp; Size</b>	<b>Size</b>	<b>Length</b>
1	SM- 389	Cardiac Pacing Wires	ELECTRODE BLUE 60CM PACING WIRE 17- 20mm 1/2 Circle Taper- point & 60mm Straight Cutting needle with Breakaway point		60cm ( ± 10cm)
2	SM- 390	Cardiac Pacing Wires	ELECTRODE BLUE 60CM PACING WIRE 26mm ½ Circle Taper- point & 90mm Straight Cutting needle with Breakaway point		60cm (±10cm)
3	SM- 395	Steel Wire	Steel Wire 1/2 Circle Cutting CCS 48 MM stainless steel Monofilament 55cm ( ± 10cm)	5	55cm (±10cm)
4	SM- 396	Steel Wire	1/2 Circle Cutting CCS 48 MM stainless steel Monofilament	6	55cm (±10cm)
5	SM- 398	Topical skin adhesive	Skin closure (2-octyl cyanoacrylate) with self-adhesive		

W

✓

			polyester mesh size 22cm, pack of 2/box, USFDA, CE approved		
6	SM- 399	Topical skin adhesive	Skin closure(2- octyl cyanoacrylate) with self adhesive polyester mesh size 60cm, pack of 2/box, USFDA, Ce approved		
7	SM- 400	Umbilical Tape	3mm x 75cm Umbilical Cotton Tape		75cm (± 10cm)
8	SM- 506	BIOPATCH	BIOPATCH Protective Disc with CHG - BIOPATCH ANTIMICROBIAL DRESSING, 1" disk (2.5 cm) w/4.0 mm center hole	2.5 cm	
9	SM- 507	BIOPATCH	BIOPATCH Protective Disc with CHG - BIOPATCH ANTIMICROBIAL DRESSING, 3/4" disk (1.9 cm) w/1.5 mm center hole	1.9 cm	
10	SM- 508	BIOPATCH	BIOPATCH Protective Disc with CHG - BIOPATCH ANTIMICROBIAL DRESSING, 1" disk (2.5 cm) w/7.0 mm center hole	2.5 cm	
11	SM- 509	INTERCEED adhesion barrier	Oxidised Regenerated Cellulose adhesion barrier - 3 X 4 " inch	3 X 4 Inch	

<b>Group- B Surgical Stapling Devices</b>			
<b>S No.</b>	<b>Item Code</b>	<b>Item</b>	<b>Tender Specs</b>
1	SM- 405	LINEAR CUTTER	Linear cutter 55/60 mm with cartridge with 4-6 rows of staple lines. Preloaded blue cartridge having closed staple height of 1.5mm. Staple row extendable approx 1.5 staples.
2	SM- 407	LINEAR CUTTER	Linear cutter cartridge Linear cutter 75/80 mm with cartridge with 4-6 rows of staple lines. With Preloaded blue cartridge having closed staple height of 1.5mm. Staple row extendable approx 1.5 staples.
3	SM- 409	LINEAR CUTTER	Linear Cutter 100 mm with cartridge with 4-6 rows of staple lines with preloaded cartridge having closed staple height of 1.5 mm
4	SM- 410 A	Linear cutter cartridge	Linear cutter cartridge 100 mm compatible with linear cutter 100 mm Blue cartridge with closed staple height of 1.5 mm
5	SM- 410 B	Linear cutter cartridge	Linear cutter cartridge 100 mm compatible with linear cutter 100 mm Green cartridge with closed staple height of 2.0mm
6	SM- 411	Endoscopic linear cutter	(normal Shaft Length) - Endoscopic linear cutter having 60 mm staple line, 57 mm cut line, with natural articulation (options of 15degree, 30 degree, 45 degree on either side), 360°rotation, 330-350 mm shaft length.
7	SM- 412	Endoscopic linear cutter	(Bariatric Shaft Length) - Endoscopic linear cutter having 60 mm staple line, 57 mm cut line, with natural articulation

✓

✓

			(options of 15°, 30°, 45° on either side), 360°rotation, 430-450mm shaft length.
8	SM-413 A	Endoscopic linear cutter cartridge	Endoscopic Linear cutter cartridge having 60 mm staple line compatible with endoscopic linear cutter 60mm, White cartridge or equivalent having closed staple height of 1.0 mm.
9	SM-413 B	Endoscopic linear cutter cartridge	Endoscopic Linear cutter cartridge having 60 mm staple line compatible with endoscopic linear cutter 60mm, blue cartridge or equivalent having closed staple height of 1.5 mm.
10	SM-413 C	Endoscopic linear cutter cartridge	Endoscopic Linear cutter cartridge having 60 mm staple line compatible with endoscopic linear cutter 60mm, Green cartridge or equivalent having closed staple height of 2.0 mm.
11	SM-415	Clip applier (preloaded) open	Multi Clip Applier (Preloaded)-Open Multi clip applier 25-30 cm with anti-clip drop ratchet mechanism, lateral & transverse grooves, 15-30° jaw angle for enhanced visibility, 20-25 preloaded titanium clip of medium size OR Multi Clip Applier (Preloaded)-Open Multiple clip applier 25cm-30cm with 15-30 preloaded titanium clips of medium size with super interlock ratchet mechanism preventing clip drop 15- 30° jaw angle for enhanced visibility.
12	SM-417	Skin stapler	Multi-directional disposable skin stapler with coated rectangular staples having dimensions of 6.5-6.9mm x 3.9-4.1mm with 35 numbers of staples with alignment indicator. Squeeze handles staple extractor
13	SM-418 A	Circular Stapler	Intra-luminal staplers with curved shaft, having controlled tissue compression/ DST through adjustable closed staple height to accommodate various tissue thicknesses 1.0-, non-grip surface, anvil detachable of cartridge sizes range from 21mm to 25 mm. OR Circular Staplers with Directional Stapling Technology with Tilt-Top Mechanism with chisel edge cutting system with safety lock-out mechanism with audible and tactile feedback and an easy to read Green "Ready to Fire" Indicator, size ranges from 21mm to 25mm.
14	SM-418 B	Circular Stapler	Intra-luminal staplers with curved shaft, having controlled tissue compression/ DST through adjustable closed staple height to accommodate various tissue thicknesses 2.0, non-grip surface, anvil detachable of cartridge sizes range from 21mm to 29 mm OR Circular Staplers with Directional Stapling Technology with Tilt-Top Mechanism with chisel edge cutting system with safety lock-out mechanism with audible and tactile feedback and an easy to read Green "Ready to Fire" Indicator, size ranges from 21mm to 28/29mm.

W

GA

15	SM-418 C	Circular Stapler	Intra-luminal staplers with curved shaft, having controlled tissue compression/ DST through adjustable closed staple height to accommodate various tissue thicknesses 2.0-, non-grip surface, anvil detachable of cartridge sizes range from 21mm to 33 mm  OR Circular Staplers with Directional Stapling Technology with Tilt-Top Mechanism with chisel edge cutting system with safety lock-out mechanism with audible and tactile feedback and an easy to read Green" Ready to Fire" Indicator, size ranges from 21mm to 33mm.
16	SM-453	Ligaclip Applier small	Clip Applicator Open Reusable Ligation Clip Applier- straight applier for small size clip 19.1cm long.
17	SM-454	LIGACLIP APPLIER 19 CM MEDIUM	Clip Applicator Open Reusable Ligation Clip Applier- straight applier for medium size clip 19.1cm long
18	SM-457	LIGACLIP ENDOSCOPIC CLIP APPLIER MEDIUM	Endoscopic Clip Applicator Endoscopic Reusable Ligation Clip Applier having alligator-jaw design, distal tip closure, rotating shaft, lateral/ transverse grooves for on vessel security, 360 rotation, for Medium size clips.
19	SM-459	LIGACLIP ENDOSCOPIC CLIP APPLIER LARGE	Endoscopic Clip Applicator Endoscopic Reusable Ligation Clip Applier having alligator-jaw design, distal tip closure, rotating shaft, lateral/ transverse grooves for on vessel security, 360 rotation, for large size clips.
20	SM-470	Endopath Reload	Endoscopic Linear cutter cartridge having 45 mm staple line compatible with Endoscopic linear cutter articulation 45 degrees ,white cartridge or equivalent having closed staple height of 1.0 mm.
21	SM-471	Endopath Reload	Endoscopic Linear cutter cartridge having 45 mm staple line compatible with Endoscopic linear cutter articulation 45 degrees , Blue cartridge or equivalent having closed staple height of 1.5 mm.
22	SM-472	Endopath Reload	Endoscopic Linear cutter cartridge having 45 mm staple line compatible with Endoscopic linear cutter articulation 45 degrees ,Gold cartridge or equivalent having closed staple height of 1.8 mm.
23	SM-473	Endopath Reload	Endoscopic Linear cutter cartridge having 45 mm staple line compatible with Endoscopic linear cutter articulation 45 degrees, Green cartridge or equivalent having closed staple height of 2.0 mm.
24	SM-474	Endopath Reload	Endoscopic Linear cutter cartridge having 45 mm staple line compatible with Endoscopic linear cutter articulation 45 degrees, Black cartridge or equivalent having closed staple height of 2.3 mm.
25	SM-478	CURVED CUTTER STAPLER 40MM	Curved Linear Cutter Curved Linear cutter having four curvilinear rows of staples in the width of 30 mm space with a 40 mm single cut line in between the staple line, simultaneous cutting and stapling, unique curved head design, lower pelvic access, manual pin placement option, reloadable cartridges

W

GA

26	SM-479	CURVED CUTTER STAPLER RELOAD	Curved Linear Cutter a) Curved blue cartridge having close staple height of 1.5mm, tactile feedback on completion of firing sequence, new anvil, knife with every cartridge
27	SM-480	CURVED CUTTER STAPLER RELOAD	Curved Linear Cutter b) Curved, green cartridge having close staple height of 2.0 mm, tactile feedback on completion of firing sequence new anvil, knife with every cartridge
28	SM-481	HAEMORRHOIDAL STAPLER PPH	Haemorrhoidal Stapler Haemorrhoidal stapler kit consists of 33mm Haemorrhoidal Circular stapler with fixed anvil
29	SM-482	POWERED LONG ARTICULATING ENDOSCOPIC LINEAR CUTTER	Endocutter Powered 60mm Articulating Endo Cutter with a disposable battery pack, for enhanced distal tip stability while firing, having closed channel in the cartridge jaw for better stability during firing, 360° rotation shaft and one handed natural articulation up to 45 degrees, Precision Machined Anvil to deliver initial, system wide compression, Wide Proximal to Distal Jaw Aperture (proximal 8mm, Distal 22mm), 3 Point Gap Control for alignment and calibration throughout the 60mm staple line, knife direction/reverse control to discontinue the firing and return the Knife, interchangeable 6 row cartridge options of closed staple height of 1.0, 1.5, 1.8, 2.0 and 2.3 mm , all fits down to 12 mm trocar sleeve, 440 mm shaft length.
30	SM-483 C	ENDOPATH STAPLER RELOADS	Endoscopic linear cutter cartridge having 60mm staple line compatible with articulating upto 45° Gold cartridge or equivalent having closed staple height of 1.8mm compatible with endoscopic linear cutter 60mm
31	SM-513		Linear Cutter 75/80 mm with cartridge with 4 rows of staple lines with mechanism for consistent staple formation, and allowing for repositioning of tissue before firing. Staple row extends approximately 1.5 staples beyond cut line. With preloaded cartridge having closed staple height of 1.5 mm OR Linear Cutter 75/80 mm with Directional Stapling Technology, integrated gap Control with 4 rows of staple lines. With preloaded cartridge having closed staple height of 1.5 mm.
32	SM-514		Linear cutter cartridge 75/80 mm compatible with above mentioned 4 rows linear cutter 75-80 mm , Cartridge with closed staple height of 1.5mm
33	SM-515		Linear cutter cartridge 75/80 mm compatible with above mentioned 4 rows linear cutter 75-80 mm , Cartridge with closed staple height of 2.0mm

W

DA

34	SM-516	<p>Linear Cutter 55/60 mm with cartridge with 4 rows of staple lines with mechanism for consistent staple formation, and allowing for repositioning of tissue before firing. Staple row extends approximately 1.5 staples beyond cut line. With preloaded cartridge having closed staple height of 1.5 mm</p> <p>OR</p> <p>Linear Cutter 55/60 mm with Directional Stapling Technology , integrated gap Control with 4 rows of staple lines. With preloaded cartridge having closed staple height of 1.5 mm.</p>
35	SM-517	Linear cutter cartridge 55/60 mm compatible with above mentioned 4 rows linear cutter 55-60 mm , Cartridge with closed staple height of 1.5mm
36	SM-518	Linear cutter cartridge 55/60 mm compatible with above mentioned 4 rows linear cutter 55-60 mm , Cartridge with closed staple height of 2.0mm
37	SM-519	<p>31 mm Circular stapler, with Rectangular Staple wire design for use even in challenging applications, detachable Slimmest Anvil Head profile of 5mm with Tilt-Top, Choice of Shaft Lengths -22 cm or 35cms, Staple Heights of 4.8 mm with closed stapled leg lengths of 1.5 mm or 2.0 mm respectively, 440 stainless steel with chisel edge cutting system with lumen/knife diameter of 22.5 mm, Slimmest Integrated trocar in shaft and an Anastomotic Lip of 1.9mm, having 30 staples, with safety lock-out mechanism with audible and tactile feedback and an easy to read Green "Ready to Fire" Indicator</p> <p>OR</p> <p>Circular Stapler Proximate Intra Luminal stapler shaft anvil detachable 31 mm should have powered firing, Atraumatic Gripping Surface Technology with 3D stapling technology</p>
38	SM-520	33 mm Circular stapler, with Rectangular Staple wire design for use even in challenging applications, detachable Slimmest Anvil Head profile of 5mm with Tilt-Top, Choice of Shaft Lengths - 22 cm or 35 cms, Staple Heights of 4.8 mm with closed stapled leg lengths of 1.5 mm or 2.0mm respectively, 440 stainless steel with chisel edge cutting system with lumen/knife diameter of 22.5mm, Slimmest Integrated trocar in shaft and an anastomotic Lip of 1.9mm, having 30 staples, with safety lock-out mechanism with audible and tactile feedback and an easy to read Green" Ready to Fire" Indicator.
39	SM-521	The Laparoscopic Gastrointestinal Stapler (Short, Slandered & Long shaft lengths), Endoscopic Stapler should have the following specifications, facility for interchanging cartridge lengths of 30 mm, 45 mm and 60mm in the same stapler. The stapler should have Tissue Gap Control Mechanism and the handle when squeezed should move only 15mm at a time, ensuring precise staple formation even at the distal end. The firing sequence of the stapler to be made possible only by activating a push button. The stapler should have a facility for 360° rotation. It should also facilitate a 22° and 45° articulation on either side of the stapler with Reticulating cartridges

40	SM-522	The Endoscopic cartridges 60MM should have the following specification: The cartridges should have 4-6 rows of staple lines irrespective of the length or color of the cartridge. Each cartridge should come with a New knife blade to help in precise cutting. Cartridges should have a wide jaw opening. (White/Blue/Green)
41	SM-523	The Endoscopic cartridges 45MM should have the following specification: The cartridges should have 4-6 rows of staple lines irrespective of the length or color of the cartridge. Each cartridge should come with a New knife blade to help in precise cutting. Cartridges should have a wide jaw opening. (White/Blue/Green)
42	SM-524	Optically guided bladeless trocar 12 mm with Optiview Technology, having wiping, wicking and absorbing technology to reduce smudging during camera reinsertion. With bilateral tissue separators, optical tip to eliminate blind entry, clear ribbed cannula, recessed stopcock valve, funnel shaped housing, duckbill secondary seal, integrated universal seal that eliminates the use of reducer in size of 100mm/150mm/180mm OR Optical bladeless trocar 12 mm with interlocking snap feature , length 100/150/180mm, optical tip to eliminate blind entry, clear universal ribbed cannula, dolphin nose tip for easier trocar placement with enhanced seal system.
43	SM-525	Optically guided bladeless trocar 5 mm with Optiview Technology, having wiping, wicking, and absorbing technology to reduce smudging during camera reinsertion. With bilateral tissue separators, optical tip to eliminate blind entry, clear ribbed cannula, recessed stopcock valve, funnel-shaped housing, duckbill secondary seal, integrated universal seal that eliminates the use of reducer in size of 100mm/150mm/180mm OR Optical bladeless trocar 5 mm with the interlocking snap feature , length 100/150/180mm, optical tip to eliminate blind entry, clear universal ribbed cannula, dolphin nose tip for easier trocar placement with enhanced seal system.
44	SM-526A	Endo Stapler 45 mm Compact Staplers battery operated and comes with a single disposable battery, Powered articulation with 55 – 60 degree articulation on either side with a total range of 110 degree articulation, One handed operation with powered continuous articulation with 3-4 buttons to articulate the jaws left, right and home position, Orientation sensor to keep articulation control consistent,, 360 degree rotation shaft, Wider jaw aperture which is 19.5mm, Precision Machined Anvil to deliver initial, system wide compression, Having closed channel in the cartridge jaw for better stability during firing, Dynamic firing for variable tissue thickness, Shaft Length of 280mm - 300mm, Compatible with 45 mm cartridge with Gripping surface Technology with options of closed staple height of 1.0, 1.5, 1.8, 2.0 and 2.3 mm, all fits

W

⊗

		down to 12mm Trocar sleeve.staple length is 45mm and cutline length 42mm
45	SM-527A	Endo Stapler 45 mm Standard Staplers battery operated and comes with a single disposable battery, Powered articulation with 55-60 degree articulation on either side with a total range of 110 degree articulation, One handed operation with powered continuous articulation with 3 buttons to articulate the jaws left, right and home position, Orientation sensor to keep articulation control consistent,, 360 degree rotation shaft, Wider jaw aperture which is 19.5mm, Precision Machined Anvil to deliver initial, system wide compression, Having closed channel in the cartridge jaw for better stability during firing, Dynamic firing for variable tissue thickness, Shaft Length 340mm-380mm, Compatible with 45 mm cartridge with Gripping surface Technology with options of closed staple height of 1.0, 1.5, 1.8, 2.0 and 2.3 mm, all fits down to 12mm-15mm trocar sleeve. Staple line length is 45 mm and cutline length is 42 mm
46	SM-528A	Endo Stapler 45 mm Long Staplers battery operated and comes with a single disposable battery, Powered articulation with 55-60 degree articulation on either side with a total range of 110 degree articulation, One handed operation with powered continuous articulation with 3-4 buttons to articulate the jaws left, right and home position, Orientation sensor to keep articulation control consistent,, 360 degree rotation shaft, Wider jaw aperture which is 19.5mm, Precision Machined Anvil to deliver initial, system wide compression, Having closed channel in the cartridge jaw for better stability during firing, Dynamic firing for variable tissue thickness, Shaft Length of 440mm-460mm, Compatible with 45 mm cartridge with Gripping surface Technology with options of closed staple height of 1.0, 1.5, 1.8, 2.0 and 2.3 mm, all fits down to 12mm-15mm trocar sleeve. Staple line length is 45 mm and cutline length is 42 mm.
47	SM-529A	Endo Stapler 60 mm Compact Staplers battery operated and comes with a single disposable battery, Powered articulation with 55-60 degree articulation on either side with a total range of 110 degree articulation, One handed operation with powered continuous articulation with 3 buttons to articulate the jaws left, right and home position, Orientation sensor to keep articulation control consistent,, 360 degree rotation shaft, Wider jaw aperture which is 22.8mm, Precision Machined Anvil to deliver initial, system wide compression, Having closed channel in the cartridge jaw for better stability during firing, Dynamic firing for variable tissue thickness, Shaft Length of 280 mm, Compatible with 60 mm cartridge with Gripping surface Technology with options of closed staple height of 1.0, 1.5, 1.8, 2.0 and 2.3 mm, all fits down to 12mm-15mm trocar sleeve. Staple line length is 60mm and cutline length is 57mm.

W

JK

48	SM-530A	Endo Stapler 60 mm Standard Staplers battery operated and comes with a single disposable battery, Powered articulation with 55 degree articulation on either side with a total range of 110 degree articulation, One handed operation with powered continuous articulation with 3 buttons to articulate the jaws left, right and home position, Orientation sensor to keep articulation control consistent,, 360 degree rotation shaft, Wider jaw aperture which is 22.8mm, Precision Machined Anvil to deliver initial, system wide compression, Having closed channel in the cartridge jaw for better stability during firing, Dynamic firing for variable tissue thickness, Shaft Length of 340 mm, Compatible with 60 mm cartridge with Gripping surface Technology with options of closed staple height of 1.0, 1.5, 1.8, 2.0 and 2.3 mm, all fits down to 12mm-15mm trocar sleeve. Staple line length is 60mm and cutline length is 57mm.
49	SM-531A	Endo Stapler 60 mm Long Staplers battery operated and comes with a single disposable battery, Powered articulation with 55 -60 degree articulation on either side with a total range of 110 degree articulation, One handed operation with powered continuous articulation with 3-4 buttons to articulate the jaws left, right and home position, Orientation sensor to keep articulation control consistent,, 360 degree rotation shaft, Wider jaw aperture which is 22.8mm, Precision Machined Anvil to deliver initial, system wide compression, Having closed channel in the cartridge jaw for better stability during firing, Dynamic firing for variable tissue thickness, Shaft Length of 440mm, Compatible with 60 mm cartridge with Gripping surface Technology with options of closed staple height of 1.0, 1.5, 1.8, 2.0 and 2.3 mm, all fits down to 12mm-15mm trocar sleeve. Staple line length is 60mm and cutline length is 57mm
50	SM 540-A	Circular Stapler battery operated having disposable battery pack with Gripping Surface technology and 3D Stapling technology with detachable anvil and adjustable closed staple height of 1.5 - 2.2 mm and open leg length of 5.2mm having staple housing outer diameter 23mm and two staple rows
51	SM 541-A	Circular Stapler Battery operated having disposable battery pack with Gripping Surface technology and 3D Stapling technology with detachable anvil and adjustable closed staple height of 1.5 - 2.2 mm and open leg length of 5.2mm having staple housing outer diameter 25mm and two staple rows
52	SM 542-A	Circular Stapler battery Operated having disposable battery pack with Gripping Surface technology and 3D Stapling technology with detachable anvil and adjustable closed staple height of 1.5 - 2.2 mm and open leg length of 5.2mm having staple housing outer diameter 29mm and two staple rows
53	SM 544-A	Linear cutter <b>75/80mm having 6 rows</b> of staple line having 6 row cartridge <b>with</b> integrated new knife blade & Linear cutter having firingknob used both sides for easy handling.
54	SM-545A	Linear cutter 55/60mm having 6 rows of staple line having 6 row cartridge with integrated new knife blade & Linear cutter having firing knob used both sides for easy handling.

W

A

55	SM-546A		Linear cutter cartridge 75/80mm having 6 rows of staple line having 6 row cartridge with integrated new knife blade & Linear cutter having firing knob used both sides for easy handling with closed staple height of 1mm 1.5mm or 2mm
56	SM-547A		Linear cutter cartridge 55/60mm having 6 rows of staple line having 6 row cartridge with integrated new knife blade & Linear cutter having firing knob used both sides for easy handling having closed staple height of 1mm 1.5mm or 2mm
57	SM-771		Reusable endoscopic powered universal stapler having powered articulation and rotation compatible with endoscopic linear cutter cartridge 30mm,45mm,60mm with complete charger and power cord.
58	SM-772		Linear Adapter Short compatible with reusable power stapling system
59	SM-773		Linear Adapter Standard compatible with reusable power stapling system
60	SM-774		Linear Adapter XL compatible with reusable power stapling system
61	SM-775		Power Control shell for reusable power stapling system

#### Group C Haemostat

Group C Haemostat				
S No	Item Code	Name of Item Haemostatic Material	Technical specification Needle Type & Size	Size
1	SM-419	Bone wax	Sterile Bone wax 2.0gm	
2	SM-420	Oxidized Regenerated Cellulose	Oxidized regenerated cellulose, thicker weave can be sutured through	3x4 inch
3	SM-421	Oxidized Regenerated Cellulose	Oxidized regenerated cellulose, thicker weave can be sutured through	6x9 inch
4	SM-422	Oxidized Regenerated Cellulose	Oxidized regenerated cellulose, based topical absorbable hemostat	4x8 inch
5	SM-423	Oxidized Regenerated Cellulose	Oxidized regenerated cellulose, based topical absorbable hemostat fibrillar/layer form fibrin material 7 layer for broad surface area coverage	2x4INCH
6	SM-424	Oxidized Regenerated Cellulose	oxidized regenerated cellulose with bacterial cidal property	2x3inch
7	SM-425	Oxidized Regenerated Cellulose	oxidized regenerated cellulose with bacterial cidal property	2x4inch
8	SM-426	Oxidized Regenerated Cellulose	oxidized regenerated cellulose	2x14inch
9	SM-427	Oxidized Regenerated Cellulose	oxidized regenerated cellulose	4x8inch
10	SM-428	Absorbable Gelatin Sponge	Absorbable Gelatin sponge should not be cross linked by formaldehyde should be sterilized by dry heat	7x5x1cm

W

S

11	SM-429	Absorbable Gelatin Matrix	Absorbable Gelatin based topical flowable hemostat with 6cc syringe prefilled with hemostatic matrix kit with thrombin
12	SM-430	Fibrin Sealant	One Vial each of BAC2 50-90mg/ml fibrinogen & Thrombin 800-1200 IU/ml human thrombin frozen solution all human formula free of bovine aprotinin triple lumen tip to prevent clogging inside the applicator maintain predictable fibrinogen viscosity in a range of temperatures. Maintain fibronectin/fibrinogen ratio with physiological ranges
13	SM-484	Intracorporeal Tissue Glues for mesh fixation	Fibrin sealant along with dual applicator for open and laparoscopic surgery both rigid 30-35cm and flexible 40-45cm

#### Group D Meshes

S No	Item Code	Name of Item	Technical specification	Size	Length
		Mesh Material	Needle Type & Size		
1	SM-431	Light weight partial absorbable mesh {Polyglecaprone + polypropylene/ Polyester + Polylactic acid (PLA)}	Light weight partial absorbable mesh {Polyglecaprone + Polypropylene/ Polyester+ Polylactic acid(PLA)}	15x15cm	± 5cm)
2	SM-432	Light weight partial absorbable mesh {Polyglecaprone + Polypropylene/Polyester + Polylactic acid(PLA)}	Light weight partial absorbable mesh {Polyglecaprone + Polypropylene/ Polyester+ Polylactic acid(PLA)}	20x15cm	± 5cm)
3	SM-433	Light weight partial absorbable mesh {Polyglecaprone + Polypropylene/Polyester + Polylactic acid(PLA)}	Light weight partial absorbable mesh {Polyglecaprone + Polypropylene/ Polyester+ Polylactic acid(PLA)}	30x15cm	± 5cm)
4	SM-435	Polypropylene/Polyester Mesh	Polypropylene/ Polyester Mesh	7.6 x 15cm	± 5cm)
5	SM-436	Polypropylene/Polyester Mesh	Polypropylene/ Polyester Mesh	15x15cm	± 5cm)
6	SM-437	Polypropylene/Polyester Mesh	Polypropylene/ Polyester Mesh	20x15cm	± 5cm)
7	SM-440	Polypropylene/Polyester Mesh	Polypropylene/ Polyester Mesh	20x20cm	± 5cm)
8	SM-441	Polypropylene/Polyester Mesh	Polypropylene/ Polyester Mesh	30x30cm	± 5cm)
9	SM-442	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/ Polyester collagen} with preplaced sutures	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/ Polyester collagen} with /Without preplaced sutures	12x12	± 5cm)

10	SM-443	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/ Polyestercoll agen} with preplaced sutures	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/Polyestercoll agen} with /Without preplaced sutures	15x10cm	± 5cm)
11	SM-444	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/Polyestercoll agen} with preplaced sutures	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/ Polyestercoll agen} with /Without preplaced sutures	15x15cm	± 5cm)
12	SM-445	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/ Polyestercoll agen} with preplaced sutures	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/ Polyestercoll agen} with /Without preplaced sutures	20x15cm	± 5cm)
13	SM-446	Tissue separating mesh {polypropylene + oxidized regenerated cellulose/ Polyestercoll agen} with preplaced sutures	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/ Polyestercoll agen} with / Without preplaced sutures	20x20cm	± 5cm)
14	SM-448	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/ Polyestercoll agen} with preplaced sutures	Tissue separating mesh {polypropylene+ oxidized regenerated cellulose/ Polyestercoll agen} with /Without preplaced sutures	30x20cm	± 5cm)
15	SM-451	Mesh fixation device	5mm single-use mesh fixation device w/25-30 absorbable Straps		
16	SM-452	Mesh fixation device	5mm single-use mesh fixation device w/ 25-30 non absorbable with titanium tacks Straps		
17	SM-621	Mesh fixation device	Powered battery operated fixation device with 30 absorbable tacks length of 7mm compatible with 5mm trocar		
18	SM-779	Anatomical Mesh	Polypropylene Anatomical Meshes with Interior flap having sealed edges around the mesh for Laparoscopic inguinal hernia repair, Right	15x10cm	± 5cm)
19	SM-780	Anatomical Mesh	Polypropylene Anatomical Meshes with Interior flap having sealed edges around the mesh for Laparoscopic inguinal hernia repair, Left	15x10cm	± 5cm)
20	SM-781	Anatomical Mesh	Polypropylene Anatomical Meshes with Interior flap having sealed edges around	16x12cm	± 5cm)

✓

✗

			the mesh for Laparoscopic inguinal hernia repair, Right		
21	SM-782	Anatomical Mesh	Polypropylene Anatomical Meshes with Interior flap having sealed edges around the mesh for Laparoscopic inguinal hernia repair, Left	16x12cm	± 5cm)
22	SM-783	Mesh fixation device	5mm absorbable, screw shaped Tackers for Mesh Fixation with 15 tacks made up of PGLA with proximal wings		

**"Bidder must quote at least 60% of the items in one group to compete for that group in the tender"**

For the ease of Bidders, fresh BOQs have been uploaded against Tender IDs mentioned against each group.

- |                      |           |                    |
|----------------------|-----------|--------------------|
| 1. Vital sutures,    | Group A1, | 2026_HME_307205_15 |
| 2. Minor Sutures,    | Group A2  | 2026_HME_307205_16 |
| 3. Miscellaneous,    | Group A3  | 2026_HME_307205_17 |
| 4. Stapling Devices, | Group-B   | 2026_HME_307205_18 |
| 5. Hemostats         | Group C   | 2026_HME_307205_19 |
| 6. Meshes            | Group D   | 2026_HME_307205_20 |

Only technical/financial bids uploaded as per revised corrigendum shall be considered for evaluation.

All the bidders are requested to keep themselves updated and submit their e- bids through e-portal [www.jktenders.gov.in](http://www.jktenders.gov.in) as per specifications & BOQs. The amendments/modifications shall be available on [www.jktenders.gov.in](http://www.jktenders.gov.in), and [www.jkmsclbusiness.com](http://www.jkmsclbusiness.com).

General Manager (K)  
JKMSCL

No. JKMSCL/Suture/Corg/686/2026/2447-51 Dated: 24-06-2026

Copy for information to the:-

1. Managing Director, JKMSCL
2. General Manager (Adm), JKMSCL.
3. Dy. General Manager(T/L) JKMSCL
4. Assistant programmer to uploaded on [www.jktenders.gov.in](http://www.jktenders.gov.in) and on [www.jkmsclbusiness.com](http://www.jkmsclbusiness.com)
5. Office Record.